	Application No.	ation No. Applicant(s)		
Notice of Allowability	09/723,823	HOLL ET AL	HOLL ET AL.	
	Examiner	Art Unit		
	Samuel P Siefke	1743		
The MAILING DATE of this communication appears on the cover sheet with the correspondence address All claims being allowable, PROSECUTION ON THE MERITS IS (OR REMAINS) CLOSED in this application. If not included herewith (or previously mailed), a Notice of Allowance (PTOL-85) or other appropriate communication will be mailed in due course. THIS NOTICE OF ALLOWABILITY IS NOT A GRANT OF PATENT RIGHTS. This application is subject to withdrawal from issue at the initiative of the Office or upon petition by the applicant. See 37 CFR 1.313 and MPEP 1308.				
1. This communication is responsive to <u>2/23/04</u> .				
2. 🔀 The allowed claim(s) is/are <u>28-30</u> .				
3. 🛮 The drawings filed on 28 November 2000 are accepted by the Examiner.				
4.				
each sheet. Replacement sheet(s) should be labeled as such in the header according to 37 CFR 1.121(d). 7. DEPOSIT OF and/or INFORMATION about the deposit of BIOLOGICAL MATERIAL must be submitted. Note the attached Examiner's comment regarding REQUIREMENT FOR THE DEPOSIT OF BIOLOGICAL MATERIAL.				
 Attachment(s) 1. ☐ Notice of References Cited (PTO-892) 2. ☐ Notice of Draftperson's Patent Drawing Review (PTO-948) 3. ☑ Information Disclosure Statements (PTO-1449 or PTO/SB/0 Paper No./Mail Date <u>5/7/01</u> 4. ☐ Examiner's Comment Regarding Requirement for Deposit of Biological Material 	6. ☑ Interview S Paper No. 8), 7. ☐ Examiner's	nformal Patent Application (PTC Summary (PTO-413), /Mail Date <u>3โนโ</u> ป) s Amendment/Comment s Statement of Reasons for Allov		

Application/Control Number: 09/723,823

Art Unit: 1743

DETAILED ACTION

Claims 28-30 are allowed.

The following is an examiner's statement of reasons for allowance: Claim 28 is allowable because the prior art does not teach or fairly suggest a sample instrument for use with a fluidic cartridge which comprises a cartridge having a first and second analysis regions; a cartridge holder for engaging said cartridge a flow cytometric measuring apparatus positioned to be optically coupled with said first analysis region said flow cytometric measuring apparatus comprising a first light source a first photodetector and a second photodetectoc and an absorption measuring apparatus positioned to be optically coupled with said second analysis region said absorption measuring apparatus comprising a second light source and a third photodetector wherein said first light source and said first photodetector are positioned to be in the same plane as a stream of particles flowing through said first analysis region.

Any comments considered necessary by applicant must be submitted no later than the payment of the issue fee and, to avoid processing delays, should preferably accompany the issue fee. Such submissions should be clearly labeled "Comments on Statement of Reasons for Allowance."

Conclusion

Any inquiry concerning this communication or earlier communications from the examiner should be directed to Samuel P Siefke whose telephone number is 571-272-1262. The examiner can normally be reached on M-F 7:00am-5:00pm.

Application/Control Number: 09/723,823 Page 3

Art Unit: 1743

If attempts to reach the examiner by telephone are unsuccessful, the examiner's supervisor, Jill A. Warden can be reached on 571-272-1700. The fax phone number for the organization where this application or proceeding is assigned is 703-872-9306.

Information regarding the status of an application may be obtained from the Patent Application Information Retrieval (PAIR) system. Status information for published applications may be obtained from either Private PAIR or Public PAIR. Status information for unpublished applications is available through Private PAIR only. For more information about the PAIR system, see http://pair-direct.uspto.gov. Should you have questions on access to the Private PAIR system, contact the Electronic Business Center (EBC) at 866-217-9197 (toll-free).

Sam P. Siefke

March 11, 2004